## MiniMuffs® Neonatal Noise Attenuators

### Instructions for Use:

Natus Medical Incorporated  
5900 First Avenue South  
Seattle  
Washington 98108 USA  
Tel.: +1-206-767-3500  
Website: [natus.com](https://www.natus.com)

EU Representative / Importer  
Natus Manufacturing Limited  
IDA Business Park  
Gort, Co. Galway, Ireland

### Associated Product Part Number:

040341

Copyright © 2021 by Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.

028196\_E 03/2021

# English

## Description:

MiniMuffs® Neonatal Noise Attenuators are single-use, soft foam earmuffs that fit around a neonate's ears to reduce the level of sound reaching the ear. The devices are sized and shaped to fit around the neonate's ear. Each pair of MiniMuffs are attached to a release liner and are removed from the liner one at a time for use on the neonates.

## Intended Use:

MiniMuffs Noise Attenuators are indicated for short-term (less than 24 hours) hospital use with neonates in a supervised setting to reduce the sound level reaching the neonate's ears.

MiniMuffs Noise Attenuators are for single use only and are disposable.

## Intended Users and Patient Target Group:

This product is for use with neonates in a hospital setting or during clinical transport to reduce noise levels and result in longer periods of sleep time.

## Clinical Benefits:

MiniMuffs Noise Attenuators reduce noise levels, which results in longer periods of sleep time, optimal oxygenation, and improved physiological stability.


## Contraindications and Side Effects:

There are no known contraindications or side effects related to Natus MiniMuffs Noise Attenuators.

## Operating Instructions:

**Important!** Follow your facility's infection control and skin care management protocols for cleaning skin, using hydrogel adhesive on premature skin, appropriate length of use, and proper disposal.


1. Remove MiniMuffs earmuffs from the liner and place them on the infant one at a time. Position in either orientation:

2. Press earmuff edges gently but firmly in place to ensure a good seal. Avoid placing earmuffs on hair. Do not press on raised center of earmuffs.

**Important!** Noise reduction is dependent on a good seal. Examine the earmuff seal when repositioning the infant.

3. To reduce adhesion of hydrogel on sensitive skin, dab hydrogel on the surface with a clean cotton swab or gauze pad prior to placing earmuff on infant.

4. Use care when removing earmuffs. Wet a cotton swab or gauze pad and slowly work under the edge of the earmuff to loosen the hydrogel. Support the underlying skin and lift the earmuff gently at the edge.

## Understanding Warning Statements:

### ⚠ WARNING

Refers to a hazardous situation that could result in death or serious injury if not avoided.

## Warnings:

### ⚠ WARNING

Be careful not to fold or bend ear. Do not force large ears into earmuffs. Do not place over compromised skin.

Examine skin according to your facility's skin care management protocol. Discontinue use if irritation, abrasion, or edema is present.

## Compliance Standards:

- ISO 10993-1: 2018 Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process

## Disposal Instructions:

Dispose of MiniMuffs Neonatal Noise Attenuators after a single use.

## Disclaimer:

Natus Medical Incorporated is not responsible for injury, infection or other damage resulting from the use of this product.

Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated and the competent authority of the Member State in which the user and/or patient is established.

Refer to the Natus website for an electronic copy of this document.

### Instructions to Access the eIFU:


A copy of the Instructions for Use in PDF format is in the associated product area:

- Newborn Care:  
<https://newborncare.natus.com/newborn-care-support>

Search for “MiniMuffs Neonatal Noise Attenuators IFU” (refer to the Product Part Number) and choose the version for your local language for the instructions to use.

The files can be printed, saved, or searched using Adobe Reader. A copy of Adobe Reader can be downloaded directly from Adobe Systems ([www.adobe.com](http://www.adobe.com)).

### Glossary of Symbols:

Symbol	Standard Reference	Standard Title	Symbol Title	Explanation
Medical Device	Not applicable	Not applicable	An indication of Medical Device	This product is a medical device.
	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices	Prescription only	Indicates that the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Not made with natural rubber latex	Indicates that the medical device is not made with natural rubber latex.
	ISO 15223-1 Symbol 5.1.1	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.
	ISO 15223-1 Symbol 5.1.3	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.1.5	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.4.3 Annex A #A.15	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Consult instructions for use	Indicates an instruction to consult an electronic instructions for use (eIFU).
	IEC 60601-1 Table D.2 #2	Medical electrical equipment – Part 1: General requirements for basic safety and Essential performance	General warning sign	Indicates a hazard of potential personal injury to patient or operator.
	MDR 2017/745	EU Medical Device Regulation	CE marking	Signifies European technical conformity.
	ISO 15223-1 Symbol 5.1.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Use-by date	Indicates the date after which the medical device is not to be used.
	ISO 15223-1 Symbol 5.2.8	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.4.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied	Do not re-use	Indicates that the medical device is intended for one use, or for use on a single patient during a single procedure.
	U.S. EPA Regulation 40 CFR Part 211, Subpart B	Code of Federal Regulations	U.S. Environmental Protection Agency logo	Products with this logo comply with the specified EPA regulation.

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